



**Product name:** smart body fat scale

**Power supply:** 3 AAA batteries

**weight range:** 5kg-180kg

**Switching:** step on scale

**weight unit:** kg

**weight accuracy:**100g

**special function:** bluetooth connect with app, low-voltage detect, overload detect

Note:Please aware that do not use with devices below, implantable medical devices, such as pacemaker and artificial lung, which are used for sustaining life, as well as wearable medical devices such as medical apparatus and ECG!

Use with above devices may cause obvious damage to body due to wrong operation and accident.




**WEIGH ONCE,  
COMPREHENSIVE ANALYSIS OF BODY CONDITION SHOWS**



weight



BMI




basic metabolism



Bone mass



rate of muscle




rate of body fat



physical age



moisture



viseral fat



protein



obesity



FFM



water content



fat weight



rate of skeletal muscle

**HOW TO USE APP**

ONLY A FEW EASY STEPS CAN CHECK ANY KINDS OF BODY INDEX

(1).Scan the QR code or search for "AIFit" to download the app.



QR CODE




AIFIT

(2).Choose your ideal unit of weight and height.



(3).Set up a new account and fill in personal information accurately. (The date of birth, your height and gender, etc.)



(4). Make sure your phone's Blue tooth® is on. Tap "My device", step on the scale with the bare feet, select the right device shown on the interface.

(5).Enter the measured interface to use. The measurement result will appear in the app interface.